

cMRI: Epilepsy Report

Epilepsy Report

Neuroradiology Group

Name Edith Epilepsy	Personal ID DEMO	Date of Birth 1958-07-01	Sex Female	Study date 2019-10-03
------------------------	---------------------	-----------------------------	---------------	--------------------------

Images

T1: Approved
QC

Volumes (Percentiles)	Right [ml]	Left [ml]
Hippocampus	3.6 (96th)	3.3 (89th)

Asymmetry (Percentiles)	Value [%]
Hippocampus	7.3 (75th)

Hippocampus - Right [ml]

Hippocampus - Asymmetry [%]

Hippocampus - Left [ml]

Comparison to cognitively normal people

Link to report: <https://usa-demo.cneuro.com/#/specialist/modules/54/cmri/65>

Report interpretation: <https://www.combinostics.com/reports/epilepsy>

Combinostics Oy, Hatanpaan valtatie 24, 33100 Tampere, Finland
<https://www.combinostics.com>

cMRI 2.1.0

1 →

2 →

3 →

4 →

← 5

← 6

← 7

← 8

- 1** Patient information and study dates.
- 7** Comparative percentile curves for Hippocampal volumes and their asymmetry.
- 2** Visual overview: T1 Segmented images.
- 8** Option to open exam in cMRI viewer.
- 3** Left and right hippocampal volumes and their percentiles and asymmetric assessment of the hippocampus and its percentile.
- 4** Link to report interpretation.
- 5** Add clinic logo or address.
- 6** Automatic quality control: T1 images.